

Sharps Disposal

Technical Information

H 270mm

W 140mm


D 145mm

Sanokil: 02 9868 4644

Sanokil's Sharps Containers are supplied in one size and collection of the units can be done on a monthly or as required basis.

Sanokil's 1.7 Litre Sharps Containers are placed in a lockable unit which enables used needles and syringes to be disposed of safely, for the protection of your customers and staff.

Sanokil Sharps Containers are an attractive design. They can be used in all public amenities or other areas as required, ideal for medical centres, hospitals, shopping centres, hotels and council public areas.



- Tamper resistant design and lockable
- Easy maintenance and servicing
- Hands free operation reduces risk of injury
- Australian standard - yellow
- Securely fitted to wall